



Ovine Plasma w/ Sodium Citrate

Collected from the source:

When searchers choose their plasma an important factor that should be taken into consideration is the source, which also emphasizes the traceability of the plasma.

Our system of vertical integration allows us to be certain of the origins and traceability of our plasma. Each manufactured batch is rigorously controlled, from the collection of plasma and throughout all stages of its treatment and production through to final packaging on our premises.

The Sheep Plasma is derived from centrifuged whole blood aseptically collected, with the addition of Sodium Citrate (8,7%).

Hemoglobin:

The hemoglobin level is measured by spectrophotometer.

pH:

The pH level is measured by phmeter.

Total Protein:

Determined by refractometer.

Storage Conditions: - 18°C

Shelf life: 4 years

Recommended use:

- Respect the product storage conditions;
- Do not use the product after the expiration date;
- Store the product in an area protected from light;
- Wear suitable clothing to handle the product to avoid contamination (e.g.: gloves, mask, hygiene cap, overall...)
- It is recommended to use the product immediately after its thaw out.